

New Product Release May 2009 - LOCATOR®



LOCATOR Light Retention Extended Range Male for divergent implants.

The addition of this new (2lbs / 907 grams) retention orange male allows you more options in providing the optimum retention choices for free-standing implant cases when divergent implants are placed into function.

New ordering information:

8915 LOCATOR Light Retention Extended Range Replacement Male (Orange) (4 Pack) \$21.00
 8915-20 LOCATOR Light Retention Extended Range Replacement Male (Orange) (20 Pack) \$95.00
 (20 Pack denotes 10% discount)

The availability of this **new light retention orange male** allows ZEST to provide you with a new LOCATOR Extended Range Male Processing Package to support your needs for divergent implants. Please consider this new 8540 package assembly which contains all 3 (green, orange, red) extended range nylon males that utilize outside retention only for divergent implants 10-20 degrees (accommodates 40 degrees between 2 or more implants).

New ordering information:







8540-2 LOCATOR Extended Range Male Processing Package (2 Pack) \$58.00
 8540-10 LOCATOR Extended Range Male Processing Package (10 Pack) \$260.00
 (10 Pack denotes 10% discount)

All LOCATOR Male Processing Packages include Titanium Denture Cap with Black Processing Male, 3 Nylon Replacement Males and Block-Out Spacer.




Retention strengths - LOCATOR Nylon Replacement Males:

Retention Values are influenced by many factors including the implant's angle; however, please use the guide below to help in determining the combined attachment retention of an overdenture.

Dual Retention 0-10 degrees	Extended Range – Angle Retention 10-20 degrees
 (8524) Clear: 5 lbs or 2268 grams	 (8547) Green: 3-4 lbs or 1361-1814 grams
 (8527) Pink: 3 lbs or 1361 grams	 (8915) Orange: 2 lbs or 907 grams
 (8529) Blue: 1.5 lbs or 680 grams	 (8548) Red: 0.5 – 1.5 lbs or 226-680 grams

Dual Retentive (inside & outside) nylon males (clear, pink, blue) should be your first choice for long-term durability and long lasting retention for minor implant angles of 0-10 degrees.



Ordering information :

8519-2 LOCATOR Male Processing Package (2 Pack) \$58.00
 8519-10 LOCATOR Male Processing Package (10 Pack) \$260.00
 (10 Pack denotes 10% discount)