

# LOCATOR® Attachment for Multi-Unit Abutments

## Flowchart for Ordering

### Free-Standing Treatment

Locator Abutment  
w/Ti Collar  
No. 8909



Abutment -Ti Alloy  
w/TiN Coating  
Collar-Ti Alloy

Processing Cap  
Male/Package  
No. 8519



Ti Cap w/Black LDPE  
Male, Block Out Spacer  
and Nylon Replacement  
Males (clear, pink, blue)

### Bar-Splinted Restoration

Locator Abutment  
w/Castable Collar  
No. 8917



Abutment -Ti Alloy  
w/TiN Coating  
Collar-Delrin

Bar Processing Cap  
Male/Package  
No. 8028



Ti Cap w/Yellow LDPE  
Male, Block Out Spacers  
and Nylon Replacement  
Males (clear, pink, blue)

CODE	LOCATOR ATTACHMENT (FREE-STANDING PLACEMENT FOR THE ANGLED POSTERIOR SITES)	PRICE
8909-2	LOCATOR Abutment w/Ti Collar for angled Multi-Unit Abutment (2 Pack)	\$250.00
8909-10	LOCATOR Abutment w/Ti Collar for angled Multi-Unit Abutment (10 Pack) 10% Discount	\$1125.00
8519-2	LOCATOR Male Processing Package (2 Pack)	\$58.00
	(Includes Denture Cap with Black Processing Male, Dual Retentive Males: Clear, Pink, Blue and Block-Out Spacer)	
8519-10	LOCATOR Male Processing Package (10 Pack) 10% Discount	\$260.00
	(Includes Denture Cap with Black Processing Male, Dual Retentive Males: Clear, Pink, Blue and Block-Out Spacer)	

CODE	LOCATOR ATTACHMENT (BAR-SPLINTED PLACEMENT)	PRICE
8917-2	LOCATOR Abutment w/Castable Collar for straight or angled Multi-Unit Abutment (2 Pack)	\$ 250.00
8917-10	LOCATOR Abutment w/Castable Collar for straight or angled Multi-Unit Abutment (10 Pack) 10% Discount	\$1125.00
8028-2	LOCATOR Yellow Bar Male Processing Package (2 Pack)	\$58.00
	(Includes Denture Cap with Yellow Processing Male, Dual Retentive Males: Clear, Pink, Blue and Block-Out Spacer)	
8028-10	LOCATOR Yellow Bar Male Processing Package (10 Pack) 10% Discount	\$260.00
	(Includes Denture Cap with Yellow Processing Male, Dual Retentive Males: Clear, Pink, Blue and Block-Out Spacer)	
	LOCATOR Implant Abutments for placement directly into the implants on the parallel anterior sites (See Table Below)	\$110.00

IMPLANT DESCRIPTION	HEIGHT OF ABUTMENT TISSUE CUFF ABOVE THE IMPLANT								
	0 mm	0.73 mm	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm
NOBELACTIVE (NP) 3.5			1401	1402	1403	1404	1405	1406	
NOBELACTIVE (RP) 4.3/5.0			1411	1412	1413	1414	1415	1416	
SELECT, REPLACE 3.5 (Purple) (Tri Lobe)	8756		8757	8758	8759	8760	8734	8735	
SELECT, REPLACE 4.3 (Gold) (Tri Lobe)	8761		8762	8763	8764	8765	8753	8754	
SELECT, REPLACE 5.0 (Blue) (Tri Lobe)			8766	8767	8768	8769	8770		
NOBEL SPEEDY REPLACE (GROOVY) 4.0	8761		8762	8763	8764	8765	8753	8754	
BRANEMARK 3.3 (NP)(HEX)		8105		8106	8107	8108	8109		
BRANEMARK 3.75, 4.0, 5.0 (RP)(HEX)		8681		8682	8683	8684	8685	8696	
BRANEMARK 5.0 (WP) (HEX)			8825 (1.1 mm)	8826	8827	8828	8829		
NOBEL SPEEDY (EXT HEX) (GROOVY) 4.0		8681		8682	8683	8684	8685	8696	
REPLACE 3.5 (HEX)			8676	8677	8678	8679	8680		
REPLACE 4.3 (HEX)			8781	8782	8783	8784	8785		
REPLACE 5.0 (HEX)			8835 (1.37 mm)	8836	8837	8838	8839		
HEX-LOCK (HL) 3.25 (HEX)			8646 (1.5 mm)	8647 (2.5 mm)	8648 (3.5 mm)	8649 (4.5 mm)	8650 (5.5 mm)		
HEX-LOCK (HL) 3.8, 4.5 (HEX)		8681		8682	8683	8684	8685	8696	
HEX-LOCK (HL) 5.0 (HEX)			8835 (1.37 mm)	8836	8837	8838	8839		
STERIOSS 3.25	(There is no abutment available for this product)								
STERIOSS 3.8	8651		8652	8653	8654	8655			

Call Zest Anchors, Inc. Customer Service at:  
**1-800-262-2310**

or place orders online at:  
**www.zestanchors.com**